


## Monoclonal anti-a1-glycoprotein From Mouse



### Method Conditions

**Column:** Cogent Bidentate C8 300™ 5µm, 300Å.  
**Catalog No.:** 40008-75P-3M  
**Dimensions:** 4.6 x 75 mm  
**Solvents:** A: DI water + 0.1% trifluoroacetic acid (TFA)  
 B: acetonitrile + 0.1% TFA  
**Gradient:**

Time (min)	%B
0.0	15
20.0	50
25.0	50
25.1	15

**Post Time:** 5 min  
**Flow Rate:** 0.5 mL/min.  
**Sample Peak:** Monoclonal anti-a1-glycoprotein, mouse 043H4848 in 0.05% w/v in DI water  
**Detection:** UV 214 nm  
**Injection:** 1 µL

### Discussion

For the last few decades glycoproteins have been a subject of interest to biochemists and biologists in many fields. These proteins are found in plasma and other biological fluids and they serve many functions in nearly every physiological process of living organisms. These macromolecules consist of a peptide chain and one or more carbohydrates linked to them. The simple gradient method used in this note for the analysis is very reproducible (%RSD about 1.5). The equilibration time between samples is very short. The peak is symmetrical and easy to integrate.

For more information visit [www.MTC-USA.com](http://www.MTC-USA.com)

Cat. No.	Description
40008-75P-3M	Cogent Bidentate C8 column for Macro Molecules, 300A, 5µm, 4.6x75mm