



Different Solvent...Different Selectivity Mechanism

Polarity or Shape Recognition – YOU CHOOSE!


Chromatogram A:



Chromatogram B:



Chromatogram C:



Steroid Compound Type

- Peak 1: Prednisolone
- Peak 2: Corticosterone
- Peak 3: Estradiol
- Peak 4: Ethinyl Estradiol
- Peak 5: Estrone
- Peak 6: Norgesterel
- Peak 7: Progesterone

Method Conditions

Mobile Phase

- Chromatogram A:** 60% Aqueous (0.1%TFA)
40% Acetonitrile (ACN)
- Chromatogram B:** 45% Aqueous (0.1%TFA) 55% MeOH
- Chromatogram C:** 56% Aqueous (0.1%TFA)
24% MeOH and 20% ACN

Column Cogent UDC-Cholesterol™

Catalog No. 69069-75R

Dimensions 75 x 4.6mm id

Flow Rate 1.0ml/minute

Detection UV, 240 nm

Temperature 15° C

Discussion of Chromatograms

An FDA requirement for Birth Control Product Analysis is the resolution of hormonal steroids ETES and ESTN.

- With Cogent Bidentate C18™ and Cogent UDC-Cholesterol™ HPLC columns, peaks ETES and ESTN are resolved with polarity-based, Reverse Phase, ACN as the organic modifier.
- With other C18 columns, these same peaks are not resolved.
- When eluting with the shape recognition in Methanol on the Cogent UDC Cholesterol™ column, ETES and ESTN are extremely well resolved. However, in this test mix ETES now co-elutes with ESTD.
- By mixing the selectivity mechanisms, polarity (ACN) and shape recognition (MeOH), partial resolution was achieved on all compounds using the UDC Cholesterol™ column. Scaling up to a 250x4.6mm column resulted in baseline resolution of all peaks.

Implications, Advantages & Uses

Using a single UDC Cholesterol™ column with orthogonal selectivity, polarity in ACN and shape recognition in Methanol can lead to unique problem solving.