

Head Space Vials & Caps

Screw Top for Better Sealing and Convenience

- Precision Threads Prevent Leaking
- Available in Different Septa Types
- Magnetic Caps w/o using Inserts
- Hand tighten with Ease, No Tools Required.



These vials and caps are a convenient and effective alternative to ordinary crimp top headspace vials and caps. Easy to use, these vials and caps are an extremely reliable system. The caps are more likely to be placed on the vial properly than crimp top vials which require crimping tools.

Built to withstand high internal pressures needed for headspace analysis, these thick wall vials are guaranteed not to burst. Since the caps are the most likely part of the system to leak, these screw caps are made with precision threads for tighter sealing and more reliability.

Save time, save operator fatigue and expect great results with these vials.

Autosampler Compatibility:

CTC Pal (Leap)	Shimadzu
Varian	Agilent
Gerstel	Others...
Atas	

For more information visit www.mtc-usa.com

Cat. No.	Description			
95010-0CV	Head Space Vials	10ml	Clear	46.0mm x 22.5mm
95010-0AV	Head Space Vials	10ml	Amber	46.0mm x 22.5mm
95020-0CV	Head Space Vials	20ml	Clear	75.5mm x 22.5mm
95020-0AV	Head Space Vials	20ml	Amber	75.5mm x 22.5mm
95020-40T-S	Cap, Butyl Red/PTFE Grey Septa		55 Shore A	1.6mm Thick.
95020-45S-T	Cap, Transparent Blue Sil/PTFE White Septa		45 Shore A	1.3mm Thick.
95020-45S-W	Cap, White Sil/PTFE Red Septa		45 Shore A	1.3mm Thick.
95020-55S-W	Cap, White Sil/PTFE Blue Septa		55 Shore A	1.6mm Thick



MTC • MicroSolv Technology Corporation

World Headquarters

1 Industrial Way West, Bldg E, Eatontown, NJ 07724 USA

Ph. 1-732-578-1777 Fx. 1-732-578-9777

www.mtc-usa.com Email: info@mtc-usa.com