

# Hot Plates, Magnetic Stirrers & Hot Plate/Magnetic Stirrers

CE



MicroSolv hotplates, stirrers and hotplate-stirrers feature an exceptionally durable, chemical resistant white ceramic work surface. Their space-efficient design (8 X 9 in. footprint), makes them ideal for use on crowded benchtops and inside of fume hoods.

Advanced microprocessor controls with convenient turn knobs allow quick, precise adjustment and maintenance of speed and temperature. Safety indicator LED's on the front panel indicate when the heating and/or stirring functions have been activated.

With a square 7.5 in work surface, all three models are compatible with a wide variety of popular sizes of borosilicate glass beakers, flasks, bottles and other vessels. The hotplate/magnetic stirrer model (55400-HS) comes complete with a support rod for mounting thermometers and temperature probes.

## Technical Data:

### Magnetic Stirrer (55400-S)

Speed Range: 60 - 1500 rpm  
 Platform: 7.5 x 7.5 in / 19 x 19 cm  
 Control: Quick Adjustment Knob  
 Dimensions: 8 (W) x 9 (D) x 4.5 (H) in.  
 Electrical Data: 120V, 60 Hz / 230V, 50/60 Hz\*

### Hotplate (55400-H)

Temperature Range: Ambient +5 to 380°C  
 Platform: 7.5 x 7.5 in / 19 x 19 cm  
 Control: Quick Adjustment Knob  
 Dimensions: 8 (W) x 9 (D) x 4.5 (H) in.  
 Electrical Data: 120V, 60 Hz / 230V, 50 / 60 Hz\*

### Hotplate/Magnetic Stirrer (55400-HS)

Speed Range: 60 - 1500 rpm  
 Temperature Range: Ambient +5 to 380°C  
 Platform: 7.5 x 7.5 in / 19 x 19 cm  
 Control: Quick Adjustment Knobs  
 Dimensions: 8 (W) x 9 (D) x 4.5 (H) in.  
 Electrical Data: 120V, 60 Hz / 230V, 50 / 60 Hz\*

## Ordering Information:

55400-S	MicroSolv Magnetic Stirrer	\$195.00
55400-S	MicroSolv Hotplate	\$195.00
55400-HS	MicroSolv Hotplate and Stirrer	\$295.00



Ph: 888-248-4972 / 732-578-1777

Fax: 732-578-9777

Web: www.mtc-usa.com