

# Autosampler Vial Caps

## The best screw cap for HPLC & LCMS



- **Clean** - Tested and Certified by LCMS
- **Clear** - New Resealing Silicone Rubber/Teflon Liner
- **Bonded** - Eliminates Contamination of the Septa by Hand, *More Secure*
- **Certified for HPLC & LCMS by LCMS**



**Clean-** The caps and septa are tested 3 times throughout the manufacturing process with very stringent pass/fail criteria that is developed for HPLC & LCMS results. Total manufacturing automation including automatically bagging the finished products in a special extractables free bag assures you that the caps are never touched by human hands that can contribute to contamination and become part of your “results”.

**Clear-** A new formulation of silicone rubber that out-performs all other septa materials for best durometer and sealability is transparent/clear. Evidence of a truly new innovation in liner technology.

**Bonded-** Only bonded caps assures you that the septa have not fallen out of the cap to be contaminated and also helps with resisting “septa plunge” when an autosampler can drive the septa into the vial.

**Pre-slit-** The caps are available in pre-slit or standard non slit to comply with your autosampler’s requirements.

**Certified-** Final LCMS QC Test Certificate of Analysis is available by lot number from our website.

Contact us for samples.



**Universal Vial Compatible-** Works with all MicroSolv Universal™ vials and other brands of 9mm, short thread vials. See our website for a complete selection of compatible vials.

Cat. No.	Description
9502S-06-LCC	LC-Clean® Cap, QC Certified SIC/PTFE, Standard,Blue, 9mm thread. 100/pack.
9502S-07-LCC	LC-Clean® Cap, QC Certified SIC/PTFE, Pre Slit Grey, 9mm thread 100/pack.



### MTC • MicroSolv Technology Corporation

#### World Headquarters

1 Industrial Way West, Bldg E, Eatontown, NJ 07724 USA

Ph. 1-732-578-1777 Fx. 1-732-578-9777

www.mtc-usa.com Email: info@mtc-usa.com

LC-Clean is a Registered Trademark of ILT

L-62

For more information visit [www.mtc-usa.com](http://www.mtc-usa.com)